# COMBILITE® SERIES

PRESSURE REGULATORS INTEGRATED  
WITH CYLINDER VALVES

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INSTRUCTION FOR USE  
BEDIENUNGSANLEITUNG  
MODE D'EMPLOI  
GEBRUIKSAANWIJZING

## 1. FOREWORD

The product complies with essential requirements of 93/42/EEC Medical Device Directive and 99/36/EEC Transportable Pressure Equipment Directive.

The valve is designed according to EN ISO 10524-3 standard.

## 2. INTENDED USE

Combilite Combination Valves are designed to be fitted to gas cylinders used for medical gases. These combination valves with a gas cylinder form medical gas packages used either as gas supply point for medical devices (anaesthetic devices, ventilating devices, incubators etc.) or for direct gas supplying to a patient's breathing mask or cannula.

GCE Combination Valves are intended to be used with the following medical gases:

- oxygen;
- nitrous oxide;
- air for breathing;
- helium;
- carbon dioxide;
- xenon;
- specified mixtures of the gases listed;
- air for driving surgical tools;
- nitrogen for driving surgical tools.

## 3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

Keep the product and its associated equipment from:

- heat sources (fire, cigarettes,, e-cigarettes,...)
- flammable materials,
- oil or grease (take a great care in the use of hand creams),
- water,
- dust.

the product and its associated equipment must be prevented from tipping over, turning over or falling.

Maintain always oxygen cleanliness standards.

Use only the product and its associated equipment in well ventilated area.

Before the first use the product shall be located in its genuine package. In case of removal from service (for transport, storage) GCE recommends using the genuine package (including inner cushioning materials).

National laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-20 / +60 °C		-30 / +60 °C
	10 / 100 %		10 / 100%
	600 / 1200 mbar		600 / 1200 mbar

**!** In case of valve storage at temperature below  $-20^{\circ}\text{C}$  do not use the valve until its temperature reaches at least  $-20^{\circ}\text{C}$ .

**!** For the valves designed to be used with mixture of gases  $\text{O}_2+\text{N}_2\text{O}$ , the lowest operating temperature is  $+5^{\circ}\text{C}$ . In common use of the valve, frosting can appear on the valve surface, which is caused by common physical reaction inside the valve when high pressure in the valve is reduced to low pressure (Joule-Thomson effect). Check that all associated equipment is connected to the valve via a hose of the length of at least 2 metres.

## 4. PERSONNEL INSTRUCTIONS

The Medical Devices Directive 93/42/EEC states that product provider must ensure that all personnel handling the product are provided with the operating instructions & performance data.

**!** Do not use the product without properly familiarization of the product and its safe operation as defined in this Instruction for use. Ensure user is aware of particular information and knowledge required for the gas in use.

## 5. PRODUCT DESCRIPTION

Combination valve combines the function of shut-off valve of a high-gas cylinder and pressure regulator for the use with medical gases. Gas from the cylinder is first controlled by the main shut-off valve and then passes through the pressure regulator and finally delivered to the patient through the flow outlet or the pressure outlet. Outlet positive pressure is fixed by the manufacturer and each valve is provided with (external or internal) low-pressure relief valve against undesirable exceeding of outlet positive pressure.

There are three basic alternatives:

- valve outlet via quick-coupler, the outlet flow is a constant one.
- valve outlet via calibrated nozzles, outlet flows can be changed by the control knob
- the combination of both alternatives selector control knob.

### A – INLET STEM

The valve is fitted to gas cylinder by a threaded inlet stem. The inlet stem can be taper threaded or parallel threaded of different size depending upon cylinder size and material.

### B – FILLING PORT


Filling port is provided for filling the gas cylinder. It includes a non-return valve (NRV). The NRV means that special filling adaptors are required to vent gas from the cylinder during the filling process.

### C – PRESSURE GAUGE

A cylinder contents pressure gauge is provided. This can be active or non-active. Versions fitted with active pressure gauge enable to check cylinder content at shut-off valve in both ON and OFF positions. The versions fitted with non-active pressure gauge enable to check cylinder content at shut-off valve in ON position only.

### D – SHUT-OFF VALVE

The combination valve is provided with a shut-off valve to isolate the gas from the cylinder from the other valve functions. It must be turned on during cylinder filling and patient therapy.

## E – RESIDUAL PRESSURE VALVE

GCE combination valve contains a residual pressure valve whose function is to retain a minimum positive pressure in the gas cylinder to avoid contamination of the cylinder content by atmospheric air.

## F, G – FLOW CONTROL HEAD “F” AND FLOW OUTLET “G” (OPTION)

GCE combination valves can be delivered with a flow control head “F”. This function is used to supply gas flow rates (l/min) at atmospheric pressure directly to a patient through the flow outlet “G”, for instance through a cannula or a facemask.

The flow outlet “G” is fitted with hose fitting (hose nipple) or a threaded type (for humidifier for instance) or with combination of both types.

## H – PRESSURE OUTLET OR QUICK-COUPLER (OPTION)

The combination valve may be fitted with a pressure outlet. The pressure outlet is the outlet direct from the low-pressure part of the valve and it is fitted with a gas specific medical quick connector also called “quick coupler” (See appendix Nr 2). The user can connect another piece of equipment to this outlet with a gas specific male probe. The quick connector self seals when the male probe is disconnected. This outlet is for supplying gas at a controlled pressure to power medical devices, for instance medical ventilator.

## EXCESS FLOW DEVICE (OPTION)

Excess flow safety valve can be provided and is located in the valve inlet stem. This provides extra safety in case of the valve breaking above the valve inlet stem (for instance at cylinder falling). It has no effect on valve function during normal service

**Note: Product colour (especially guard, flow control knob and shut-off valve) might not follow any gas colour coding.**

# 6. USE OF PRODUCT

## 6.1. DISTRIBUTOR OPERATIONS

### 6.1.1. BEFORE USE

The hereunder operations are performed for instance by the provider (operator) of the device, trained hospital staff or emergency care staff or trained supplier ensuring delivery of complete medical package to the user, patient.

#### VISUAL INSPECTION BEFORE USE

- Check if there is visible damage on the medical gas primary package (including product labels and marking). If it shows signs of external damages, remove the product from service and suitably identify its status.
- Visually check if the medical gas primary package is contaminated; if needed, apply the cleaning procedure hereunder.
- Ensure turning-off at the shut-off valve by turning the hand wheel in a clockwise direction to ZERO position. Do not use excessive force (Maximum recommended closing torque is 2 Nm)
- Set flow control knob on ZERO position - ensure that the flow control knob engages correctly.

Remove caps from flow outlet and/or pressure outlet. Keep caps in a safe place for reuse during transport or storage.

#### LEAK-TIGHT CHECK BEFORE USE

- Set flow control knob (if any) on ZERO position - ensure the flow control knob engages correctly.
- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.

- **Visually and by sound check possible leakages from:**

- Flow outlet
- Filling port
- Pressure gauge
- Shut-off valve
- The lower part of flow control knob (fuse venting)

 **If any leakage is detected, apply the procedure described in Chapter 7c. and return the valve to GCE for service.**

#### FUNCTIONAL CHECK BEFORE USE

- Set flow control knob on ZERO position - ensure that flow control knob engages correctly.
- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

 **Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Check if pressure is indicated on pressure gauge. If pressure indicator enters red zone, return the medical package (cylinder with valve) for its refilling.
- Check that there is a gas flow at each set positions in both clockwise and anticlockwise turning direction (for instance by sound or checking presence of bubbles in a humidifier).
- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force (Maximum recommended closing torque is 3 Nm).
- Reset flow control knob on ZERO position when it's not venting anymore - until the flow control knob engages.
- For combination valve fitted with pressure outlet, ensure it is in working condition by connecting and disconnecting quick-coupler probe (refer to the section hereunder for pressure outlet connection).

### 6.1.2. USER OUTLET(S) CONNECTION

#### LIST OF RECOGNIZED ACCESSORIES

##### *To be connected to the flow outlet:*

- Humidifier extension (the outlet hose nipple is in horizontal position, a "L-shaped" extension may be necessary to use a humidifier, ask for special extension & humidifiers from GCE),
- Oxygen tubing connected with breathing mask or cannula,
- Gas saver,
- Atomizer.

##### *To be connected to the pressure outlet:*

- Low-pressure hoses,
- Flow meters,
- Venturi suction ejectors.

 **Before connecting any accessory or medical device to the regulator, always check that they are fully compatible with the product connection features & the product performances.**

#### PRESSURE OUTLET

- Ensure that the quick-coupler probe is compatible with the pressure outlet features.
- Connect the quick-coupler probe.
- Ensure that quick-coupler probe is well engaged.

 **If the pressure indicator enters red zone and the gas is supplied via pressure outlet, outlet pressure P2 can be lower than 3.6 bar (see appendix nr. 1 for each product).**

 **If pressure outlet is to be connected to medical device that requires high gas flow (for instance pulmonary ventilator that requires gas flow 100 l/min at the minimum pressure 2.8 bar), check required flow of connecting medical device with pressure and flow characteristics of the valve stated in appendix nr. 1. For assurance of sufficient performance (pressure and flow characteristics of the valve) do not use medical package if pressure indicator enters red zone.**

## FLOW OUTLET

**Before connecting any accessory to the flow outlet make sure that the patient is not connected before using the valve.**

- Ensure that the hose/humidifier is compatible with the flow outlet features.
- Push hose onto the regulator flow outlet/screw-on the humidifier.
- Ensure that the hose/humidifier is well engaged/screwed-on.

### 6.1.3. AFTER USE

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force.
- Ensure the flow control knob (if any) is on ZERO position.
- Disconnect all downstream equipment from user outlets.
- Refit pressure outlet and flow outlet protection cap. Before refitting the caps, ensure they are clean.

## 6.2. OPERATIONS PERFORMED BY THE USER / PATIENT

### 6.2.1. USE OF USER OUTLET CONNECTIONS

#### USE OF PRODUCT THROUGH THE FLOW OUTLET (FLOW SETTING)

- Ensure that the flow control knob is on ZERO position
- Ensure that the accessory is connected to the flow outlet.
- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

**Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Set flow control knob on one required of the available flow rates.

**Always ensure that the flow control knob has been engaged and not placed between two settings otherwise the flow control head will not deliver correct or any flow of medical gas.**

**The same variants of flow control head can have an “end stop” in between the maximum flow position and the zero position. Do not try to apply an excessive torque on the flow control knob when it stops on the maximum flow position (during clockwise rotation) or in zero position (anti-clockwise rotation).**

**Medical gas flow rate must be prescribed by a doctor.**

#### AFTER COMPLETION OF THE THERAPY

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force.
- Vent gas pressure from downstream equipment.
- Reset flow control knob on ZERO position when it's not venting anymore - until the flow control knob engages
- Disconnect the tube from the flow outlet.

#### USE OF PRODUCT THROUGH PRESSURE OUTLET

- Ensure that the flow control knob is on ZERO position (if any).
- Ensure that the accessory IS NOT connected to the pressure outlet.
- Open slowly the combination shut-off valve by turning the hand wheel in anticlockwise direction approx 1 to 1½ turns.

**Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.**

- Connect the accessory to the pressure outlet.
- After completion of the therapy
- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force.
- Vent gas pressure from downstream equipment.
- Disconnect the accessory from the pressure outlet.

## 6.2.2. AFTER USE

- Turn off at the shut-off valve by turning the hand wheel in a clockwise direction to stop position. Do not use excessive force.
- Ensure that the flow control knob is on ZERO position (if any).
- Disconnect all downstream equipment from user outlets.
- Refit pressure outlet and flow outlet protection cap. Before refitting the caps, ensure they are clean.

## 7. ACCESSORY

### BED HANGER

- Bed hanger is used with a guard. It is designed for suspension of medical package – cylinder and combination valve with the guard – onto the patient's bed. The maximum hanger load is with the use of a 5-litre cylinder.

### GCE FILLING ADAPTORS

**GCE recommends using GCE filling adaptors; Non GCE filling adaptors can damage the non-return valve in the filling port and impact the safety of the product.**

## 8. CLEANING

Remove dirt with a soft cloth damped in oil free soap water & rinse with clean water. Used cleaning agents shall comply with oxygen cleanliness standards.

Disinfection can be carried out with an alcohol-based solution (spray or wipes).

If other cleaning solutions are used, check that they are not abrasive and they are compatible with the product materials (including labels) and gas (convenient cleaning solution - i.e. Meliseptol).

**Do not use cleaning solutions containing ammonia!**

**Do not expose to water or any other liquid.**

**Do not expose to high temperature (such as autoclave).**

**To apply the cleaning solution do not spray it as the spray may enter into the inner parts of combination valve and cause contamination or damage.**

**Do not use pressure wash as it could damage or contaminate the combination valve.**

**If the inner parts of the combination valve have been contaminated do not continue to use the combination valve under any circumstances. It must be withdrawn from service.**

## 9. SERVICE, PRODUCT LIFE TIME AND MAINTENANCE

### 9.1. SERVICE PERIOD AND PRODUCT LIFE TIME

#### SERIAL NUMBER AND MANUFACTURE DATE

A nine-figure serial number stamped on valve body consists of the data as follows:

YY MM XXXXX

YY: year of manufacture

MM: month of manufacture

XXXXX: sequence number of product

For example: Serial number 050300521 indicates the valve manufactured in 2005, in March, with sequence number 521.

An eleven-figure serial number stamped on valve body consists of the data as follows:

YYYY MM XXXXX

YYYY: year of manufacture

MM: month of manufacture

XXXXX: sequence number of product

For example: Serial number 20050300521 indicates the valve manufactured in 2005, in March, with sequence number 521.

## PRODUCT LIFE TIME AND WASTE MANAGEMENT

Maximum life time of the product is 10 years from the manufacturing date.

At the end of the product's life time (10 years maximum), the product must be withdrawn from service. The provider of the device shall prevent the reuse of the product and handle the product in compliance with "Directive of European Parliament and Council 2008/98/EC on waste".

## 9.2. MAINTENANCE, REPAIR AND SERVICE OPERATIONS

### PRE/POST FILLING CHECKS

Perform functional and leak-tight checks before and after each filling. See Annex 3 for further instructions. It should not be assumed that Pre / Post Filling Checks (refer to Appendix nr 3) cover every safety procedure or method required by local regulations, or that abnormal or unusual circumstances may not warrant or suggest further requirements or additional procedures. Contact GCE if further information is needed.

### REPAIR OPERATIONS

Repairs consist of corrective operations (i.e.: involve product known as faulty).

Repair activities of GCE combination valves must be carried out by GCE or authorized repair centres. Contact GCE

for further information about service available in your area.

Repair activities can generally be performed with the combination valve still fitted on the cylinder.

Repair activities for all GCE combination valve families, which do not require an authorization for repair, are the replacements of the following damaged components:

- guard,
- labels,
- caps and detachable hose nipples.

Depending on the combination valve type, there are also other replacements/repairs of components, which do not require repair authorisation. Contact GCE for further information.

 **All labels on the equipment must be kept in good, legible condition by the provider and the user during the entire product life time.**

 **Use only genuine GCE components.**

Any product sent back to GCE (or GCE authorised centre) for repair or maintenance shall properly be packaged according GCE recommendations. Contact GCE for further information about available maintenance packaging set.

The purpose of the maintenance shall clearly be specified (repair, overall maintenance). For product to be repaired, the fault short description and any reference to a claim nr should be indicated.

## 10. GLOSSARY

	Consult Instruction for use.		Suitable for Hospital care use
	Caution!		Suitable for Home care use
	Keep away from heat and flammable material		Suitable for Emergency care use
	Keep away from oil and grease		Serial number

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Tschechische Republik      © GCE s.r.o